



material : black polyoxymethylene (POM)
weight : 0,85 kg

TARGETS

- allows to shut down a bag ring or glove ring on confinement containers from class 1 to 4

FEATURES

- can be assembled on any type of ferrule of internal diameter 186 mm
- two flat sides allow introduction of the door inside the glove box

ASSEMBLING

- possibility of assembling from inside or outside the container
- eccentric clamping device

CONSUMABLES

- o-ring

CERTIFICATE

- certificate of conformity supplied with the part (in the packaging)

Port	ID	OD	L	Part ID
Ø156mm	156mm	220mm	68mm	GBA99000278
Ø186mm	186mm	237mm	67mm	GBA99000279
Ø220mm	217mm	270mm	70mm	GBA99000280
Ø250mm	249mm	207x312mm	68mm	GBA99000281